

FIRECR+ MEDICAL READERS



FireCR+ Medical Reader Product Line

FireCR+ Medical Readers are becoming the go-to products for radiology departments and specialty practices looking for high-quality, versatile, highly reliable CR solutions.

FireCR+ Medical-80 Reader – Up to 70 plates/hour (IP 35cm x 43cm)

FireCR+ Medical-60 Reader – Up to 60 plates/hour (IP 35cm x 43cm)

FireCR+ Medical-40 Reader – Up to 40 plates/hour (IP 35cm x 43cm)

FireCR+ Medical-20 Reader – Up to 20 plates/hour (IP 35cm x 43cm)

- Operates two universal cassettes with imaging plate sizes of 24cm x 30cm & 35cm x 43cm
- Select from four models depending on your imaging speed requirements
- Upgrade as your imaging volume increases
- Compact & lightweight—choose tabletop or space-saving wall-mount option for space-challenged practices
- Transport case on wheels offers portability
- Features QuantorMed+ Imaging Software—next generation workstation software with an easy-to-use interface for intuitive & fast operation



Vertical Wall Mount

- Provides a footprint of less than 5 inches
- Allows installation in clinics with limited space

FireCR+M Medical Reader Product Line

New to the 3DISC Imaging product family, the FireCR+M Medical Reader is capable of scanning multiple cassette sizes—including the most-used imaging plates. The tabletop system offers the same high image quality, speed, and upgradeability that are hallmarks of 3DISC's CR reader product line.

FireCR+M Medical-80 Reader – Up to 70 plates/hour (IP 35cm x 43cm)

FireCR+M Medical-60 Reader – Up to 60 plates/hour (IP 35cm x 43cm)

FireCR+M Medical-40 Reader – Up to 40 plates/hour (IP 35cm x 43cm)

FireCR+M Medical-20 Reader – Up to 20 plates/hour (IP 35cm x 43cm)

- Scans a full range of the most-used cassette sizes of 18 cm x 24 cm, 24 cm x 30 cm & 35 cm x 43 cm
- The FireCR+M Medical Reader's 24 cm x 30 cm size cassettes are perfect for podiatry examinations
- Select from four models depending on your imaging speed requirements
- Upgrade as your imaging volume increases
- Compact, tabletop system
- Features QuantorMed+ Imaging Software—next generation software with an easy-to-use interface for intuitive & fast operation



Transport case

- All FireCR+ and FireCR+ M Readers are delivered in a solid transport case on wheels, providing clinics with a portable solution